



FEATRUES

- Linux-based system
- Supporting dual display
- Connectivity of upto 32 bed-side patient monitors
- Wired(Ethernet or serial port) and wireless connection
- Bi-direction communication
- Mass storage of patient data

Color Central Station System DS 8000

FEATURE



▲ Major Monitoring



▲ Main Screen



▲ Dynamic Waveform Recall

SPECIFICATION

CONNECTIVITY :

- 1) Wired serial(base) : Serial connection
- 2) Wired LAN(Optional) : RJ45 LAN connection
- 3) Wireless WLAN(Optional) : 801.11.b.g wireless access point connection

DISPLAY :

- Max 16beds x 2monitor
- Waveform : ECG, SpO2, Resp., IBP, EtCO2
- Anesthetic Agent Data : HR, PR, ST, RR, NIBP, Temp, SpO2, IBP, EtCO2, CO
- Special bed monitoring window
- Drug dose calculation

ALARM :

- Audible and visual
- Upper / Low alarm limits
- Arrhythmia alarm

- 96 hours real-time recall, multi-channel full disclosure recall
- 500 alarm events auto storage and recall
- 1000 Arrhythmia events storage and recall
- 1000 sets measurement tubular data
- 10000 patient data storage



▲ Trend graphic

SPECIFICATIONS SUBJECT TO BE CHANGED WITHOUT PRIOR NOTICE



20610 MANHATTAN PL. #112
TORRANCE, CA 90501, USA

TEL +1-310-541-8421, 324-1350

FAX +1-310-943-2271

E-mail sales@dst-inc.com

URL www.dst-inc.com